


PRODUCT SHEET

Date of issue: 31/01/23

Disinfection mop for Sanety tool, mixed cotton



REF.	NOTES	NET WEIGHT (kg)	SIZE (cm) ↓ ↓ ↓	PCK	VOLUME (m³)	GROSS WEIGHT (kg)
8136		0,08	65x15	100	0,055	8,73

APPLICATIONS

Disinfection mop head to be used with Sanety tool.

ECO BOX



ECO-FRIENDLY PRODUCT

Product sold in packaging compliant with Legislative Decree 116 of 03-09-2020 concerning the obligation of environmental labeling for correct disposal at the end of use. 99% of the cardboard packaging used by the company is made on average with 80% recycled paper

MATERIAL

50% cotton - 50% polyester (+/- 10%)

LIFE CYCLE

Washable max 200 cycles if washing instructions and warnings are respected

CARE AND MAINTENANCE

WARNINGS

Follow carefully the washing instructions and warnings on the detergent packaging (temperature and dosing) before washing
 We recommend the use of normal detergent for textiles, no fabric softener
 Always a pre-wash cycle in cold water without detergent to remove any residue of chemical solution
 We suggest to use the specific net bags for the mop washing

WASHING INSTRUCTIONS



MAX washing temperature: 95° C



Bleaching possible - Do not use chlorine or bleach



Drying-chamber



No ironing

Disclaimer: Filmop International Srl reserves the right to update the contents of this document at any time and without notice.



PRODUCT SHEET

Date of issue: 31/01/23

⊗ Do not dry clean

Do not use softener

Use alkaline detergents with Ph < 11

RELATED PRODUCTS



SANETY - disinfection tool with bottle 1,5 L

GENERAL PRODUCT SAFETY

The products sold by us do not need a product certification, as they are in compliance with the definition in Legislative Decree 206 of 06.09.2005

DECLARATION OF CONFORMITY

All the products sold by us are in compliance with the laws in force regarding safety and are covered by an insurance policy "Civil responsibility for products"